

Histamine Conjugate

Catalog No.	AG005	Quantity:	1 mg
Description:	Histamine Conjugate: Histamine-G-Pc (Protein Carrier). G = Glutaraldehyde.		
Applications:	Immunohistochemistry, immunocytochemistry.		
Storage & Handling:	Lyophilized antigen are stable at least 1 year. For reconstitution, use sterile ddH ₂ O. This solution is stable at +4°C for 15 days.		
Related Product(s):	Catalog No. PA2023 : Anti-conjugated Histamine (rabbit). Catalog No. PA2024 : Anti-conjugated Histamine (rat).		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

