

L-Tyrosine Conjugate

Catalog No.	AG008	Quantity:	1 mg
Description:	L-Tyrosine Conjugate: L-Tyrosine-G-Pc (Protein Carrier). G = Glutaraldehyde.		
Applications:	Immunohistochemistry, immunocytochemistry.		
Storage & Handling:	Lyophilized antigen are stable at least 1 year. For reconstitution, use sterile ddH ₂ O. This solution is stable at +4°C for 15 days.		
Related Product(s):	Catalog No. PA2028 : Anti-conjugated Tyrosione (rabbit). Catalog No. PA2029 : Anti-conjugated Tyrosine (rat).		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

