

L-Phenylalanine Conjugate

| | | | |
|--------------------------------|---|------------------|------|
| Catalog No. | AG009 | Quantity: | 1 mg |
| Description: | L-Phenylalanine Conjugate: L-Phenylalanine-G-Pc (Protein Carrier). G = Glutaraldehyde. | | |
| Applications: | Immunohistochemistry, immunocytochemistry. | | |
| Storage & Handling: | Lyophilized antigen are stable at least 1 year. For reconstitution, use sterile ddH ₂ O. This solution is stable at +4°C for 15 days. | | |
| Related Product(s): | Catalog No. PA2030 : Anti-conjugated Phenylalanine (rabbit). | | |

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

