

## L-glutamine Conjugate

|                                |   |                  |      |
|--------------------------------|---|------------------|------|
| <b>Catalog No.</b>             | AG017   | <b>Quantity:</b> | 1 mg |
| <b>Description:</b>            | L-glutamine Conjugate: L-glutamine-G-Pc (Protein Carrier).<br>G = Glutaraldehyde.   |                  |      |
| <b>Applications:</b>           | Immunohistochemistry, immunocytochemistry.  |                  |      |
| <b>Storage &amp; Handling:</b> | Lyophilized antigen are stable at least 1 year.<br>For reconstitution, use sterile ddH <sub>2</sub> O. This solution is stable at +4°C for 15 days. |                  |      |
| <b>Related Product(s):</b>     | Catalog No. <a href="#">PA2039</a> : Anti-conjugated L-glutamine (rat).   |                  |      |

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

