

L-iso-leucine Conjugate

Catalog No.	AG019	Quantity:	1 mg
Description:	L-iso-leucine Conjugate: L-iso-leucine-G-Pc (Protein Carrier). G = Glutaraldehyde.		
Applications:	Immunohistochemistry, immunocytochemistry.		
Storage & Handling:	Lyophilized antigen are stable at least 1 year. For reconstitution, use sterile ddH ₂ O. This solution is stable at +4°C for 15 days.		
Related Product(s):	Catalog No. PA2041 : Anti-conjugated L-iso-leucine (rabbit).		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

