

Trans-hydroxyproline Conjugate

Catalog No.	AG022	Quantity:	1 mg
Description:	Trans-hydroxyproline Conjugate: Trans-hydroxyproline-G-Pc (Protein Carrier). G = Glutaraldehyde.		
Applications:	Immunohistochemistry, immunocytochemistry.		
Storage & Handling:	Lyophilized antigen are stable at least 1 year. For reconstitution, use sterile ddH ₂ O. This solution is stable at +4°C for 15 days.		
Related Product(s):	Catalog No. PA2044 : Anti-conjugated Trans-hydroxyproline (rabbit).		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

