

D-Cysteine Conjugate

Catalog No. AG027 **Quantity:** 1.0 mg

Description: This product is the D-Cysteine conjugated to glutaraldehyde and BSA

Formulation: Lyophilized.

Reconstitution: **Centrifuge vial prior to opening.** Reconstitute product to desired concentration in sterile water containing 0.1% merthiolate.

Storage & Stability: Lyophilized antigen is stable for at least two years. Reconstituted solution is stable at +4°C for one month in darkness.

Applications: ELISA, Immunohistochemistry, immunocytochemistry.

The optimal concentration should be determined by the user for each specific application.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

