

D-Tryptophan Conjugate

Catalog No. AG029 **Quantity:** 1.0 mg

Description: This product is the D-Tryptophan conjugated to glutaraldehyde and BSA

Formulation: Lyophilized.

Reconstitution: **Centrifuge vial prior to opening.** Reconstitute product to desired concentration in sterile water containing 0.1% merthiolate.

Storage & Stability: Store lyophilized product at 4°C; this product is stable for 2 years. The reconstituted product is stable at 4°C for 2 months.

Applications: Immunohistochemistry, immunocytochemistry.

The optimal concentration should be determined by the user for each specific application.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

