

## L-dihydroxyphenylalanine (L-DOPA) Conjugate

<b>Catalog No.</b>	AG032	<b>Quantity:</b>	1 mg
<b>Description:</b>	L-DOPA Conjugate: L-DOPA-G-Pc (Protein Carrier). G = Glutaraldehyde.		
<b>Applications:</b>	Immunohistochemistry, immunocytochemistry.		
<b>Storage &amp; Handling:</b>	Lyophilized antigen are stable at least 1 year. For reconstitution, use sterile ddH <sub>2</sub> O. This solution is stable at +4°C for 15 days.		
<b>Related Product(s):</b>	Catalog No. <a href="#">PA2003</a> : Anti-conjugated L-dihydroxyphenylalanine (L-DOPA) (rabbit). Catalog No. <a href="#">MA2001</a> : Anti-conjugated L-dihydroxyphenylalanine (L-DOPA) (mouse).		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

