

5-Hydroxytryptamine Conjugate

Catalog No.	AG044	Quantity:	1 mg
Description:	5-Hydroxytryptamine conjugate: 5-hydroxytryptamine-G-Pc (Protein Carrier). G = Glutaraldehyde		
Formulation:	Lyophilized antigen		
Applications:	Immunohistochemistry, immunocytochemistry.		
Reconstitution:	Centrifuge vial prior to opening. For reconstitution, use distilled H ₂ O.		
Storage & Stability:	Lyophilized antigen is stable at least 1 year. Reconstituted solution is stable at +4°C for 15 days.		
Related Product(s):	Catalog No. PA2012 : Anti-conjugated 5-hydroxytryptamine (rabbit)		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences, Inc.
480 Neponset Street
Bldg 12A
Canton, MA 02021

Toll Free: 888-769-1246
Phone: 781-828-0610
Fax: 781-828-0542

E-mail: info@cellsciences.com
Web site: www.cellsciences.com