

Tryptamine Conjugate

Catalog No. AG045 **Quantity:** 1.0 mg

Description: Tryptamine conjugate: Tryptamine-G-BSA (G = Glutaraldehyde; BSA = Bovine Serum albumin)

Formulation: Lyophilized

Reconstitution: **Centrifuge vial prior to opening.** Reconstitute with distilled water containing 0.1% merthiolate.

Storage & Stability: Lyophilized antigen are stable at least 2 year at -20°C. When reconstituted according to the above directions, product is stable at 2-4°C for 2 months. **Avoid repeated freeze-thaw cycles.**

Applications: Immunohistochemistry, immunocytochemistry.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

