

## 5-Methoxytryptamine Conjugate

**Catalog No.** AG046 **Quantity:** 1.0 mg

**Description:** This product is the 5-Methoxytryptamine conjugated to glutaraldehyde and BSA

**Formulation:** Lyophilized.

**Reconstitution:** **Centrifuge vial prior to opening.** Reconstitute product to desired concentration in sterile water containing 0.1% merthiolate.

**Storage & Stability:** Lyophilized antigen is stable at least 2 years. Reconstituted solution is stable at +4°C for 1 month.

**Applications:** Immunohistochemistry, immunocytochemistry, ELISA

The optimal concentration should be determined by the user for each specific application.

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

