

4,5-Dihydroxytryptamine Conjugate

Catalog No. AG049 **Quantity:** 1.0 mg

Description: This product is the 4,5-Dihydroxytryptamine conjugated to Glutaraldehyde and BSA

Formulation: Lyophilized.

Reconstitution: **Centrifuge vial prior to opening.** Reconstitute product to desired concentration in sterile water containing 0.1% merthiolate.

Storage & Stability: Lyophilized antigen is stable at least 2 years. Reconstituted solution is stable at +4°C for 1 month.

Applications: Immunohistochemistry, immunocytochemistry, ELISA

The optimal concentration should be determined by the user for each specific application.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.