

Malondialdehyde-PC Conjugate

Catalog No.	AG057	Quantity:	1.0 mg
Description:	Malondialdehyde conjugated to carrier protein.		
Formulation:	Lyophilized		
Reconstitution:	For reconstitution, use pure water with 0.1% merthiolate. For better solubility, 10 µL of 1N sodium hydroxide may be necessary.		
Applications:	Immunohistochemistry, immunocytochemistry.		
Storage & Stability:	Upon receipt, store desiccated at -20 °C for up to 1 year. Store reconstituted product at 4 °C for up to 1 month in darkness. Product is light-sensitive.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

