

Malondialdehyde Conjugate

Catalog No.	AG057	Quantity:	1.0 mg
Description:	This product is the Malondialdehyde conjugated to Bovine Serum Albumin (BSA)		
Formulation:	Lyophilized		
Reconstitution:	For reconstitution, use pure water with 0.1% merthiolate.		
Applications:	Immunohistochemistry, immunocytochemistry. The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Lyophilized antigen is stable at least 2 years. The reconstituted product is stable at 4°C for 1 month in darkness.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.