

## 5,7-Dihydroxytryptamine Conjugate

**Catalog No.** AG068 **Quantity:** 1.0 mg

**Description:** This product is the 5,7-Dihydroxytryptamine conjugated to Glutaraldehyde and Bovine Serum Albumin (BSA)

**Formulation:** Lyophilized.

**Reconstitution:** **Centrifuge vial prior to opening.** Reconstitute product to desired concentration in sterile water containing 0.1% merthiolate.

**Storage & Stability:** Lyophilized antigen is stable at least 2 years. The reconstituted product is stable at 4°C for 1 month in darkness.

**Applications:** Immunohistochemistry, immunocytochemistry, ELISA

The optimal concentration should be determined by the user for each specific application.

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

