

Nitroso-L-DOPA Conjugate

Catalog No. AG080 **Quantity:** 1.0 mg

Description: This product is the Nitroso-L-DOPA conjugated to Nitrosylated, Glutaraldehyde, and Bovine Serum Albumin (BSA)

Formulation: Lyophilized.

Reconstitution: **Centrifuge vial prior to opening.** Reconstitute product to desired concentration in sterile water containing 0.1% merthiolate.

Storage & Stability: Lyophilized antigen is stable at least 2 years. The reconstituted product is stable at 4°C for 1 month in darkness.

Applications: Immunohistochemistry, immunocytochemistry
The optimal concentration should be determined by the user for each specific application.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

