

## Nitroso-Tyramine Conjugate

**Catalog No.** AG081 **Quantity:** 1.0 mg

**Description:** This product is the Nitroso-Tyramine conjugated to Nitrosylated, Glutaraldehyde, and Bovine Serum Albumin (BSA)

**Formulation:** Lyophilized.

**Reconstitution:** **Centrifuge vial prior to opening.** Reconstitute product to desired concentration in sterile water containing 0.1% merthiolate.

**Storage & Stability:** Lyophilized antigen is stable at least 2 years. The reconstituted product is stable at 4°C for 1 month in darkness.

**Applications:** Immunohistochemistry, immunocytochemistry  
The optimal concentration should be determined by the user for each specific application.

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**