

Agmatine ((4-Aminobutyl)guanidine) Conjugate

Catalog No.	AG085	Quantity:	1.0 mg
Description:	This product is the Agmatine conjugated to Glutaraldehyde, and Bovine Serum Albumin (BSA)		
Formulation:	Lyophilized.		
Reconstitution:	Centrifuge vial prior to opening. Reconstitute product to desired concentration in sterile water containing 0.1% merthiolate.		
Storage & Stability:	Lyophilized antigen is stable at least 2 years. The reconstituted product is stable at 4°C for 1 month in darkness.		
Applications:	Immunohistochemistry, immunocytochemistry The optimal concentration should be determined by the user for each specific application.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

