

## Auxin (Indole-3-Acetic Acid) Conjugate

<b>Catalog No.</b>	AG086	<b>Quantity:</b>	1.0 mg
<b>Description:</b>	This product is the Auxin (Indole-3-Acetic Acid) conjugated to Bovine Serum Albumin (BSA)		
<b>Formulation:</b>	Lyophilized.		
<b>Reconstitution:</b>	<b>Centrifuge vial prior to opening.</b> Reconstitute product to desired concentration in sterile water containing 0.1% merthiolate.		
<b>Storage &amp; Stability:</b>	Lyophilized antigen is stable at least 2 years. The reconstituted product is stable at 4°C for 1 month in darkness.		
<b>Applications:</b>	Immunohistochemistry, immunocytochemistry The optimal concentration should be determined by the user for each specific application.		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

