

Abscissic Acid Conjugate

Catalog No. AG088 **Quantity:** 1.0 mg

Description: Abscissic acid conjugate: abscissic acid-BSA (BSA = Bovine Serum albumin)

Formulation: Lyophilized

Reconstitution: **Centrifuge vial prior to opening.** Reconstitute with distilled water containing 0.1% merthiolate.

Applications: Immunohistochemistry, immunocytochemistry.
The optimal concentration should be determined by the user for each specific application.

Storage & Stability: Lyophilized antigen is stable at least 2 years. The reconstituted product is stable at 4°C for 1 month in darkness.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

