

## Glutathione-Protein carrier (Pc) Antigen

**Catalog No.** AG094 **Quantity:** 1.0 mg

**Description:** Glutathione conjugate: Glutathione-G-Pc. G = Glutaraldehyde

**Reconstitution:** **Centrifuge vial prior to opening.** Reconstitute in pure water to the desired concentration.

**Storage & Stability:** Lyophilized antigen is stable at least 1 year at -20°C. Antigen reconstituted in distilled water is stable at 4°C for 15 days. **Avoid freeze-thaw cycles.**

**Applications:** Immunocytochemistry

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**