

## alpha-Tocopherol (Vitamin E) Conjugate

|                                 |   |                  |        |
|---------------------------------|---|------------------|--------|
| <b>Catalog No.</b>              | AG100   | <b>Quantity:</b> | 1.0 mg |
| <b>Description:</b>             | This product is the alpha-Tocopherol (Vitamin E) conjugated to Bovine Serum Albumin (BSA)   |                  |        |
| <b>Formulation:</b>             | Lyophilized.  |                  |        |
| <b>Reconstitution:</b>          | <b>Centrifuge vial prior to opening.</b> Reconstitute product to desired concentration in sterile water containing 0.1% merthiolate.    |                  |        |
| <b>Storage &amp; Stability:</b> | Lyophilized antigen is stable at least 2 years. The reconstituted product is stable at 4°C for 1 month in darkness.                     |                  |        |
| <b>Applications:</b>            | Immunohistochemistry, immunocytochemistry.<br>The optimal concentration should be determined by the user for each specific application. |                  |        |

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

