

Vitamin H/D-biotin Conjugate

Catalog No.	AG101	Quantity:	1.0 mg
Description:	Vitamin H/D-biotin Conjugate: Vitamin H/D-biotin-BSA (BSA = Bovine Serum albumin)		
Formulation:	Lyophilized		
Reconstitution:	Centrifuge vial prior to opening. Reconstitute with distilled water containing 0.1% merthiolate.		
Storage & Stability:	Lyophilized antigen are stable at least 2 year at -20°C. When reconstituted according to the above directions, product is stable at 2-4°C for 2 months. Avoid repeated freeze-thaw cycles.		
Applications:	Immunohistochemistry, immunocytochemistry.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

