

Nicotinamide Conjugate

Catalog No.	AG105	Quantity:	1.0 mg
Description:	Niacin (also known as Vitamin B3, Nicotinic Acid and Vitamin PP). The liver can synthesize Niacin from the essential amino acid Tryptophan. Nicotinamide is conjugate to Bovine Serum Albumin.		
Formulation:	Lyophilized		
Reconstitution:	Centrifuge vial prior to opening. Reconstitute with distilled water containing 0.1% merthiolate.		
Applications:	Immunohistochemistry Immunocytochemistry The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Lyophilized antigen is stable for 2 years at -20°C. Reconstituted antigen is stable for 2 months at 2-4°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

