

## Mouse IgG1 Isotype Control (Clone B-Z1) Azide Free mAb

<b>Catalog No.</b>	CDM377	<b>Quantity:</b>	200 µg
<b>Description:</b>	The mouse monoclonal antibody is for use as an isotype control.		
<b>Concentration:</b>	1 mg/ml		
<b>Conjugate:</b>	Unconjugated		
<b>Specificity:</b>	Not specific for any known human antigens.		
<b>Host:</b>	Mouse		
<b>Isotype:</b>	IgG1		
<b>Immunogen:</b>	Rat monocytes and macrophages.		
<b>Clone:</b>	B-Z1		
<b>Formulation:</b>	Liquid in PBS. Sterile-filtered and treated to remove endotoxins. Carrier and preservative free.		
<b>Applications:</b>	Flow Cytometry, Functional Assays		
<b>Application Notes:</b>	Negative control for biological tests. The optimal concentration should be determined by the user for each specific application.		
<b>Storage &amp; Stability:</b>	Stable at 2-8°C for 1 year. For long term storage, freeze in aliquots at -20 to -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

