

Mouse IgG1 Isotype Control (Clone B-Z1) PE mAb

Catalog No.	CDM379	Quantity:	100 tests
Description:	The mouse monoclonal antibody is for use as an isotype control.		
Conjugate:	PE		
Specificity:	Not specific for any known human antigens.		
Host:	Mouse		
Isotype:	IgG1		
Immunogen:	Rat monocytes and macrophages.		
Clone:	B-Z1		
Formulation:	Lyophilized from PBS + 1% BSA + 0.1% sodium azide. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Reconstitution:	Reconstitute with 1 ml sterile deionized water.		
Applications:	Flow Cytometry		
Application Notes:	Use 10 μ l to label 10^6 cells or 100 μ l of whole blood. The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Stable at 2-8°C for 6 months after reconstitution. DO NOT FREEZE.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®

65 Parker Street

Unit 11

Newburyport, MA 01950

Toll Free: 888-769-1246

Phone: 978-572-1070

Fax: 978-992-0298

E-mail: info@cellsciences.com

Website: www.cellsciences.com