

Mouse IgG2a Isotype Control (Clone B-Z2) Azide Free mAb

Catalog No.	CDM383	Quantity:	200 µg
Description:	The mouse monoclonal antibody is for use as an isotype control.		
Concentration:	1 mg/ml		
Conjugate:	Unconjugated		
Specificity:	Not specific for any known human antigens.		
Host:	Mouse		
Isotype:	IgG2a		
Immunogen:	Rat monocytes and macrophages.		
Clone:	B-Z2		
Formulation:	Liquid in PBS. Sterile-filtered and treated to remove endotoxins. Carrier and preservative free.		
Applications:	Flow Cytometry, Functional Assays		
Application Notes:	Negative control for biological tests. The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Stable at 2-8°C for 1 year. For long term storage, freeze in aliquots at -20 to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.