


Rat Anti-Mouse IgG_{2b} Clone SB74g PE mAb

Catalog No.	CMM037	Quantity:	200 µg
Description:	Rat Anti-Mouse IgG _{2b} Clone SB74g PE mAb		
Concentration:	100 µg/mL		
Specificity:	Mouse IgG _{2b}		
Immunogen:	Mouse IgG _{2b} hybridoma		
Isotype:	Rat (Sprague Dawley) IgG _{2b} λ		
Clone:	SB74g		
Conjugate:	PE (R-phycoerythrin)		
Formulation:	Liquid volume of 1.0 mL in a buffer of PBS/NaN ₃ and a stabilizing agent.		
Isotype Control	Rat IgG _{2b} -PE		
Applications:	Applications for relevant formats of this clone include FLISA and Flow Cytometry.		
Application Notes:	Working dilutions; FLISA ≤ 1 µg/mL. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C. DO NOT FREEZE		

FLISA plate was coated with purified mouse IgG1, IgG2a, IgG2b, IgG3, IgM, and IgA. Immunoglobulins were detected with serially diluted Rat Anti-Mouse IgG2b-PE.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

