

F12

Rabbit Anti-Mouse Factor XII IgG Fraction HRP pAb

Catalog No.	CPM007A CPM007B	Quantity:	1 mg 10 mg
Alternate Names:	HAF, Hageman factor, FXII		
Description:	HRP conjugated Rabbit polyclonal IgG fraction to mouse factor XII.		
Gene ID:	58992		
Specificity:	Recognizes mouse Factor XII		
Host:	Rabbit		
Immunogen:	Recombinant full length wild-type protein from insect cells (glycosylated)		
Isotype:	IgG		
Conjugate:	HRP		
Formulation:	Frozen liquid in 0.05 M PBS + 1% BSA + 0.01% Thimerosal		
Purification:	IgG fraction purified by immobilized protein A.		
Applications:	ELISA, WB		
Application Notes:	The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	When stored at -80°C, product is stable for 5 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

