

F12

Rabbit Anti-Mouse Factor XII Affinity Purified Biotin pAb

Catalog No.	CPM009A CPM009B	Quantity:	100 µg 1.0 mg
Alternate Names:	HAF, Hageman factor, FXII		
Description:	Biotin conjugated Rabbit affinity purified polyclonal IgG fraction to mouse factor XII.		
Gene ID:	58992		
Specificity:	Recognizes mouse Factor XII		
Host:	Rabbit		
Immunogen:	Recombinant full length wild-type protein from insect cells (glycosylated)		
Isotype:	IgG		
Conjugate:	Biotin		
Formulation:	Frozen liquid in 0.05M Sodium Phosphate + 0.1 M NaCl + 1 mM EDTA, pH 6.6.		
Purification:	Affinity purified by immobilized mouse factor XII		
Applications:	ELISA, WB		
Application Notes:	The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	When stored at -80°C, product is stable for 5 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

