

Corn Trypsin Inhibitor

Catalog No.	CRT186B CRT186C	Quantity:	1 mg 10 mg
Description:	Corn Trypsin Inhibitor		
Concentration:	4.3 mg/ml		
Source:	Fresh sweet corn		
Molecular Weight:	12.5 kDa		
Formulation:	Liquid in 0.02 M Tris-HCl + 0.15 M NaCl, pH 7.4		
Storage & Stability:	Store in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

