

Cyno Monkey Plasma Albumin

Catalog No.	CS162A	Quantity:	10 mg
	CS162B		50 mg

Description: Albumin is a water-soluble protein with considerable structural stability which makes up 60 percent of the total protein of plasma. It functions as a carrier of hormones, enzymes, fatty acids, metal ions, and medicinal products. Prepared from cynomolgus monkey plasma by proprietary fractionation methods and lyophilized. >95percent pure by SDS-PAGE. Add deionized water or buffer to original or desired volume, aliquot and freeze unused portion.

Concentration: Lyophilized

Source: Cyno monkey plasma

Molecular Weight: 66 kDa


Formulation: Salt free lyophilized powder

Purity: >95% by SDS-PAGE analysis

Storage & Stability: When stored at -80°C, product is stable for three years from date of delivery.

Volume: 1.0 ml before lyophilization

1. Cyno Monkey Albumin (10 µg) NR
2. Molecular Weight Markers



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.