

## Bovine IgG Protein G Purified Isotype Control

<b>Catalog No.</b>	CS208A CS208B	<b>Quantity:</b>	50 mg 100 mg
<b>Description:</b>	Purified from normal serum by immobilized Protein G. >95 percent pure by SDS-PAGE and preservative free.		
<b>Concentration:</b>	5.0 mg/ml		
<b>Purity:</b>	IgG fraction		
<b>Specificity:</b>	This antibody has been selected as an isotype control as it has no known specificity.		
<b>Host:</b>	Bovine		
<b>Isotype:</b>	IgG		
<b>Conjugate:</b>	Protein G Purified		
<b>Formulation:</b>	Liquid in 0.05M Sodium Phosphate + 0.1M NaCl + 1mM EDTA, pH 6.6		
<b>Purification:</b>	>95 percent pure by SDS-PAGE and preservative free.		
<b>Applications:</b>	ELISA, WB		
<b>Storage &amp; Stability:</b>	When stored at -80°C, product is stable for 5 years from date of delivery. <b>Avoid repeated freeze-thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

