

PLG Dog Plasminogen

Catalog No.	CS211A CS211B	Quantity:	1.0 mg 10 mg
Description:	Prepared from fresh US origin dog (canine) plasma by immobilized lysine chromatography		
Concentration:	10.0 mg/ml		
Gene ID:	403602		
Source:	Dog plasma		
Molecular Weight:	93 kDa		
Formulation:	Liquid in 0.1M HEPES; 0.1M NaCl; pH 7.4		
Purity:	>95% by SDS-PAGE analysis		
Storage & Stability:	When stored at -80°C, product is stable for three years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

