

## ALB

### Horse Plasma Albumin

<b>Catalog No.</b>	CS214A	<b>Quantity:</b>	100 mg
	CS214B		200 mg

**Description:** Albumin is a water-soluble protein with considerable structural stability which makes up 60 percent of the total protein of plasma. It functions as a carrier of hormones, enzymes, fatty acids, metal ions, and medicinal products. Prepared from US origin equine (horse) plasma by proprietary fractionation methods and lyophilized. >95 percent pure by SDS-PAGE. Add deionized water or buffer to original or desired volume, aliquot and freeze unused portion.

**Concentration:** Lyophilized

**Gene ID:** 100034206

**Source:** Horse plasma

**Molecular Weight:** 66 kDa

**Formulation:** Lyophilized powder Salt Free

**Purity:** >95% by SDS-PAGE analysis

**Storage & Stability:** When stored at -80°C, product is stable for three years from date of delivery. **Avoid repeated freeze-thaw cycles.**

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

