

Rat Liver Lyophilized Powder

Catalog No.	CS243A CS243B	Quantity:	1 g 2 g
Description:	Sprague Dawley rat livers are frozen, lyophilized, ground and bottled.		
Concentration:	Lyophilized		
Source:	Rat liver		
Applications:	The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	When stored at -80°C, product is stable for 5 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.