

PLAT tPA Coated Plate

Catalog No.	CS406A	Quantity:	5 plates
	CS406B		25 plates

Description: 96 well plate with 8 removable strips coated with recombinant human tissue-type Plasminogen Activator (tPA) at 100 ul/well and blocked with 300 ul/well. Ideal for irreversible binding of active Plasminogen Activator Inhibitor (PAI-1) for functional ELISA experiments. Human tPA forms a covalent complex with PAI-1 of all species including human, mouse, rat, porcine and rabbit.

Concentration: 0.005 mg/ml

Source: CHO cells

Form: Plate

Storage & Stability: When stored at 2-4°C, product is stable for 12 months from date of delivery.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

