

## Armenian Hamster Immunoglobulin G Isotype Control Clone HTK888 Biotin

<b>Catalog No.</b>	CSI11925 CSI11926	<b>Quantity:</b>	50 µg 200 µg
<b>Description:</b>	This antibody was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse, rat and human tissues.		
<b>Concentration:</b>	0.5 mg/ml		
<b>Immunogen:</b>	Trinitrophenol + KLH		
<b>Isotype:</b>	Armenian Hamster IgG		
<b>Clone:</b>	HTK888		
<b>Formulation:</b>	Phosphate-buffered solution + pH 7.2, containing 0.09% sodium azide. <b>Precaution:</b> Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
<b>Purification:</b>	The antibody was purified by affinity chromatography, and conjugated with biotin under optimal conditions. The solution is free of unconjugated biotin.		
<b>Applications:</b>	FC, ICFC		
<b>Recommended Usage:</b>	Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis as negative control. Use at concentrations comparable to those of the specific antibody of interest.		
<b>Storage &amp; Stability:</b>	The antibody solution should be stored undiluted at 4 °C. <b>Do not freeze.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

