

## TNF


### Armenian Hamster Anti-Mouse TNF-alpha Clone TN3-19.12 LE/NA mAb

<b>Catalog No.</b>	CSI12347 CSI12348	<b>Quantity:</b>	50 µg 0.5 mg
<b>Alternate Names:</b>	Tumor necrosis factor-α, Cachectin, Necrosin, Macrophage cytotoxic factor (MCF), Differentiation inducing factor (DIF), TNFSF-2, TNF-alpha,		
<b>Description:</b>	TNF-α is secreted by macrophages, monocytes, neutrophils, T-cells (principally CD4+), and NK-cells. Many transformed cell lines also secrete TNF-α. Monomeric mouse TNF-α is 156 amino acid protein (N-glycosylated) with a reported molecular weight of 17.5 kD protein. TNF-α forms multimeric complexes; stable trimers are most common in solution. A 26 kD membrane form of TNF-α has also been described. TNF-α binding to surface receptors elicits a wide array of biologic activities including: cytolysis and cytostasis of many tumor cell lines in vitro, hemorrhagic necrosis of tumors in vivo, increased fibroblast proliferation, and enhanced chemotaxis and phagocytosis in neutrophils. The TN3-19.12 antibody reacts with mouse, rat, and rabbit tumor-necrosis factor-α (TNF-α) proteins. The TN3-19.12 antibody can neutralize the bioactivity of natural or recombinant TNF-α.		
<b>Concentration:</b>	1.0 mg/ml		
<b>Gene ID:</b>	21926		
<b>Structure:</b>	TNF superfamily; dimer/trimer; 18-150 kD (Mammalian).		
<b>Regulation:</b>	Type II integral membrane protein processed by TACE for secretion; upregulated by interferons, IL-2, GM-CSF, substance P, bradykinin, PAF, immune complexes, cyclooxygenase; downregulated by IL-6, TGF-β, vitamin D3, prostaglandin E2, PAF antagonists.		
<b>Host:</b>	Armenian Hamster		
<b>Immunogen:</b>	<i>E. coli</i> -expressed, recombinant mouse TNF-α		
<b>Isotype:</b>	Armenian Hamster IgG		
<b>Clone:</b>	TN3-19.12		
<b>Bioactivity:</b>	Paracrine/endocrine mediator of inflammatory and immune functions; selectively cytotoxic for transformed cells; endothelial cell alterations; chemoattractant.		
<b>Formulation:</b>	Phosphate-buffered solution, pH 7.2, containing no preservative. 0.2 µm filter sterilized. Endotoxin level is < 0.1 EU/µg of the protein (< 0.01 ng/µg of the protein) as determined by the LAL test.		
<b>Purification:</b>			



The LE/NA (Low Endotoxin, Azide-Free) antibody was Purified by affinity chromatography.

- Receptors:** TNFRSF1A (TNF-R1, CD120a, TNFR-p60 Type  $\beta$ , p55); TNFRSF1B (TNF-R2, CD120b, TNFR-p80 Type A, p75)
- Reactivity:** Mouse, Rat, Cross-Reactivity: Rabbit (Lapine)
- Applications:** ELISA Capture, ELISPOT Capture, ICFC, IP, Neut, WB
- Recommended Usage:** Each lot of this antibody is quality control tested by ELISA assay. For ELISPOT applications, a concentration range of 4-8  $\mu\text{g/ml}$  is recommended. For ELISA capture applications, a concentration range of 2-6  $\mu\text{g/ml}$  is recommended. To obtain a linear standard curve, serial dilutions of mouse TNF- $\alpha$  recombinant protein ranging from 500 to 4  $\text{pg/ml}$  are recommended for each ELISA plate. It is recommended that the reagent be titrated for optimal performance for each application.
- Storage & Stability:** The antibody solution should be stored undiluted at 4 °C. This LE/NA solution contains no preservative; handle under aseptic conditions.
- Cellular Sources:** Activated monocytes, neutrophils, macrophages, T cells, B cells, NK cells, LAK cells
- Cellular Targets:** Monocytes, neutrophils, macrophages, T cells, fibroblasts, endothelial cells, osteoclasts, adipocytes, astroglia, microglia.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

