

Recombinant Bovine RabGDI

Catalog No.	CSI14870	Quantity:	50 µg
Description:	Rab GDP Dissociation Inhibitor (RabGDI) forms complexes with cytosolic prenylated GDP-bound Rab proteins and delivers them to the target membrane. It also retrieves GDP-bound Rab proteins from the membrane releases.		
Source:	Sf9 insect cells		
Formulation:	Liquid. Supplied in 40 mM Hepes, pH 7.2 + 50 mM NaCl + 100 µM GDP and 3 mM DTE.		
Purity:	≥90% by SDS-PAGE		
Storage & Stability:	Stable for 1 year in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

