

SERPINF2

A2AP, Rabbit Anti-Human AntiPlasmin IgG Fraction pAb

Catalog No.	CSI19726A	Quantity:	1.0 mg
	CSI19726B		10 mg

Description: Alpha 2-AntiPlasmin (A2AP) is a serine protease inhibitor (serpin) responsible for inactivating plasmin, an important enzyme that participates in fibrinolysis and degradation of various other proteins. This protein is encoded by the SERPINF2 gene.

Concentration: 10.0 mg/mL

Volume: 0.1 mL

Molecular Weight: 160 kDa

Host: Rabbit Serum

Formulation: Frozen Liquid in 0.05 M Sodium Phosphate pH 6.6 + 0.1 M NaCl + 1 mM EDTA

Purification: Protein A purified

Applications: ELISA, WB.
The optimal concentration should be determined by the user for each specific application.

Storage & Stability: Store at -70°C.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

