

Ecotin (E. coli serine protease inhibitor)

Catalog No.	CSI20044A CSI20044B	Quantity:	1.0 mg 5 mg
Description:	Serine protease inhibitor purified from E. coli. At least 98percent pure as determined by SDS-PAGE.		
Concentration:	2.5 mg/ml		
Source:	E. coli		
Molecular Weight:	16,096 (monomer)		
Formulation:	Liquid in 1mM HCl		
Purity:	98% by SDS-PAGE analysis		
Storage & Stability:	When stored at -80°C product is stable for 3 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.