

## Phe-Pro-Arg-CMK Fluorescein labeled

<b>Catalog No.</b>	CSI20045A CSI20045B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Description:</b>	Phe-Pro-Arg-CMK Fluorescein labeled		
<b>Concentration:</b>	10 mg/ml		
<b>Source:</b>	Chemical		
<b>Molecular Weight:</b>	848 Da		
<b>Formulation:</b>	Frozen liquid containing DMSO		
<b>Purity:</b>	>95% by HPLC analysis		
<b>Endotoxin Level:</b>	< 0.1 ng/µg of protein		
<b>Storage &amp; Stability:</b>	Store at -80°C. Product is stable for 3 years from delivery date. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

