

Glu-Gly-Arg-Chloromethylketone

Catalog No.	CSI20047A CSI20047B	Quantity:	5.0 mg 25.0 mg
Description:	Glu-Gly-Arg-Chloromethylketone		
Source:	Chemical		
Molecular Weight:	466 Da		
Formulation:	Lyophilized from a solution containing 1 mM HCl.		
Purity:	>95% by HPLC analysis		
Endotoxin Level:	< 0.1 ng/μg of protein		
Specific Activity:	1 mole of Glu-Gly-Arg-Chloromethylketone inhibits 0.34 moles Human Factor Xa		
Storage & Stability:	Store at -80°C. Product is stable for 3 years from delivery date. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

