

Phe-Phe-Arg-Chloromethylketone

Catalog No.	CSI20061A	Quantity:	10 mg
	CSI20061B		100 mg

Description: Phe-Phe-Arg-CMK (FFRCK). Potent and specific irreversible inhibitor of plasma and glandular kallikreins. Inhibits the cleavage of atrial natriuretic peptide (ANP 1-126) by kallikreins. Also inhibits adrenal kininogenase, a kallikrein-like enzyme.

Source: Chemical

Molecular Weight: 501 Da

Formulation: Lyophilized powder, salt free

Purity: >95% by HPLC analysis

Endotoxin Level: < 0.1 ng/μg of protein

Storage & Stability: Store at -80°C. Stable for 3 years from delivery. **Avoid repeated freeze-thaw cycles.**

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.