

Mouse Anti-Bovine Serum Albumin Clone 9E2C2 mAb

Catalog No.	CSI20162A CSI20162B	Quantity:	1 mg 10 mg
Description:	Mouse Anti-Bovine Serum Albumin Clone 9E2C2 mAb		
Concentration:	10.0 mg/ml		
Host:	Mouse		
Immunogen:	Full length native, purified bovine serum albumin		
Isotype:	IgG1 kappa		
Clone:	9E2C2		
Formulation:	Frozen liquid buffered in 0.05 M Sodium Phosphate + 0.1 M NaCl + 1 mM EDTA, pH 6.6.		
Purification:	IgG fraction purified by immobilized Protein A.		
Applications:	ELISA, Western Blot		
Storage & Stability:	Store at -80°C.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.