

Donkey IgG, Protein G Purified

Catalog No.	CSI20166A CSI20166B	Quantity:	50 mg 100 mg
Description:	Donkey IgG purified by immobilized Protein G. IgG fraction, purified from normal serum by immobilized Protein G.		
Concentration:	10 mg/ml		
Source:	Donkey serum		
Molecular Weight:	160 kDa		
Formulation:	Frozen liquid buffered in 0.02 M Sodium Phosphate + 0.15 M NaCl, pH 7.4. Preservative free.		
Purity:	>95% pure by SDS-PAGE.		
Endotoxin Level:	< 0.1 ng/μg of protein		
Specific Activity:	This antibody has been selected as an isotype control as it has no known specificity.		
Storage & Stability:	Store at -80°C.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

