

Glu-Gly-Arg-Chloromethylketone Biotin

Catalog No.	CSI20170A CSI20170B	Quantity:	1.0 mg 5.0 mg
Description:	Biotinylated Glu-Gly-Arg-CMK is a biotin labeled peptide to be used as a specific, irreversible inhibitor of Factor Xa.		
Concentration:	10 mg/ml		
Source:	Chemical synthesis		
Molecular Weight:	882 Da		
Formulation:	Frozen liquid in 0.01 M HCl solution		
Purity:	>95% by HPLC analysis		
Endotoxin Level:	< 0.1 ng/μg of protein		
Specific Activity:	1 mole of PPACK inhibits 1 mole Human Thrombin		
Storage & Stability:	Store at -80°C. Stable for 3 years from delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

