

## SERPINE1

### Recombinant Mouse PAI-1 Latent Fraction Biotin

<b>Catalog No.</b>	CSI20199A CSI20199B	<b>Quantity:</b>	0.5 mg 1.0 mg
<b>Alternate Names:</b>	PAI1, Planh1, Serpine1, serine (or cysteine) peptidase inhibitor		
<b>Description:</b>	The latent fraction of mouse recombinant PAI-1 is biotin labeled on lysine residues.		
<b>Concentration:</b>	2.5 mg/ml		
<b>Gene ID:</b>	18787		
<b>Source:</b>	<i>E. coli</i>		
<b>Molecular Weight:</b>	43 kDa		
<b>Formulation:</b>	Frozen liquid in 0.1 M Tris-HCl Buffer + 0.15 M NaCl, pH 7.4		
<b>Purity:</b>	>95% by SDS-PAGE analysis		
<b>Endotoxin Level:</b>	< 0.1 ng/μg of protein		
<b>Storage &amp; Stability:</b>	Store at -80°C. Stable at -80°C for 3 years from delivery. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

